Retinoblastoma Detector ©

![Python](https://img.shields.io/badge/Python-3.10-blue)

![PyTorch](https://img.shields.io/badge/Framework-PyTorch-red)

![Streamlit](https://img.shields.io/badge/UI-Streamlit-green)![License:

MIT](https://img.shields.io/badge/License-MIT-yellow)

An AI-powered application designed to detect Retinoblastoma, a rare eye cancer, from images of eyes or faces. The project uses deep learning, computer vision, and a user-friendly web app interface for diagnosis and research purposes.



- 1. **Eye Detection**: Automatically detects eyes in face photos using Mediapipe.
- 2. **Disease Classification**: Classifies eyes as Healthy or Retinoblastoma and provides confidence scores.
- 3. **Interactive Web App**: Built with Streamlit for user-friendly interaction.
- 4. **Customizable Model**: Trainable Convolutional Neural Network (CNN).

Setup Instructions

1. Clone the Repository

Clone the repository to your local machine:

```bash

git clone https://github.com/your-username/retinoblastoma-detector.git cd retinoblastoma-detector

#### 2. Install Dependencies

This project uses Poetry for dependency management. Install Poetry if you don't already have it:

"bash

pip install poetry

..

Activate the environment and install all dependencies:

```bash

poetry install

•

3. Prepare the Training Data

Organize the `Training/` folder with the following structure:

٠.,

Ensure that there are sufficient images in both categories to train the model effectively.

4. Train the Model

Run the training script:

"bash

poetry run python train.py

This script:

- 1. Loads the training data from the `Training/` directory.
- 2. Trains the Convolutional Neural Network (CNN).
- 3. Saves the trained model to the 'models/' directory as 'retinoblastoma_detector.pth'.

(#) Usage Instructions

- 1. **Upload an Image**: Choose an image of an eye or face. Supported formats: `.jpg`, `.jpeg`, `.png`.
- 2. **View Predictions**: For face images, eyes are detected and cropped automatically. Each eye is classified as Healthy or Retinoblastoma, and confidence scores are displayed.
- 3. **Adjust Confidence Threshold**: Use the sidebar slider to filter predictions by confidence levels.

How It Works

- 1. **Eye Detection**: Uses Mediapipe to detect and crop eyes from face images.
- 2. **Image Preprocessing**: Resizes images to `224x224` and normalizes pixel values using mean `[0.5]` and standard deviation `[0.5]`.
- 3. **Model Architecture**: A CNN is used for binary classification (Healthy vs. Retinoblastoma).
- 4. **Inference**: Cropped eye regions are passed through the trained model, and predictions are output with confidence scores.

Dependencies

All dependencies are listed in 'requirements.txt'. Notable ones include:

- **torch**: Deep learning framework.
- **torchvision**: For data augmentation and image preprocessing.
- **opency-python-headless**: Image processing library.
- **mediapipe**: For eye detection in face images.
- **streamlit**: For building the interactive web app.

Future Improvements

- 1. **Pretrained Models**: Use transfer learning with pretrained architectures like ResNet.
- 2. **Dataset Expansion**: Include more diverse images for better generalization.
- 3. **Batch Processing**: Support uploading and processing multiple images simultaneously.
- 4. **Mobile Optimization**: Enhance the Streamlit interface for better mobile usability.
- 5. **API Integration**: Add RESTful API endpoints for remote inference.

Contributing

We welcome contributions to enhance the project! To contribute:

- 1. **Fork the repository**.
- 2. Create a new branch:

```
```bash
```

git checkout -b feature/your-feature

3. Commit your changes:

```bash

git commit -m "Add your feature"

4. Push your changes:

```bash

git push origin feature/your-feature

5. Open a \*\*Pull Request\*\* on GitHub.



#### **License**

This project is licensed under the \*\*MIT License\*\*. See the `LICENSE` file for details.

## **Contact**

For questions or feedback, feel free to reach out:

- \*\*Your Name\*\*
- [Your GitHub Profile](https://github.com/your-username)
- [Your Email Address] (mailto:your-email@example.com)